


**U-CyTech B.V.**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
+31 85 073 1460  
info@ucytech.com  
www.ucytech.com

## Data sheet      Human IFN- $\gamma$ /IL-10 FluoroSpot kit

**Catalogue no.**      CT850-F2 (2-plate format)

**Species reactivity**      human



*Typical example of IFN- $\gamma$ /IL-10 specific spots produced by  $2.5 \times 10^4$  human PBMC*  
In green: IFN- $\gamma$ , in red: IL-10, in yellow: IFN- $\gamma$ /IL-10.

Preincubation:      42h  
Incubation FluoroSpot:      20h  
Stimulation:      PMA/ionomycin

*References can be found in the Reference Database at [www.ucytech.com/references](http://www.ucytech.com/references).*

