


U-CyTech B.V.
Yalelaan 48
3584 CM Utrecht
The Netherlands
+31 85 073 1460
info@ucytech.com
www.ucytech.com

Data sheet Human IgM B cell ELISPOT kit

Catalogue no. CT782-PR5 (5-plate format)

Species reactivity human



Total number of IgM producing cells
1x10⁴ cells/well
Stimuli preincubation: IL-2 + R848

Example of IgM specific spots produced by human PBMC

B cell activation: 4 days preincubation (2x10⁶ PBMC/ml) in presence of recombinant IL-2
and R848.

Incubation ELISPOT: 24 hours

References can be found in the Reference Database at www.ucytech.com/references.

